Insulin Cooler

7.78 OZ **EN**

PRODUCT FEATURES

Reliable Cooling Tube

Stable Vacuum Insulation

Compact and Portable

Smaller Tubes
Larger Capacity

Liner 304 Stainless Steel

15°C (59°F) Gel Phase
Change Tube

HOW TO USE

- **Cool the Cooler Cup:** Refrigerate at ≤40°F (4°C) for 3 hours or freeze at 0°F (-18°C) for 1 hour.
- **Cool the Ice Tubes:** Refrigerate at ≤40°F (4°C) for 90 minutes or freeze at 0°F (-18°C) for 45 minutes.
- **Insert Insulin:** Place the insulin pen/vial in the pre-cooled cup, add ice tubes, and secure the lid tightly.
- **Transport:** Avoid shaking or impacts. Open the lid only when needed, and promptly put the insulin pen back to cup after use.
- **Ice Tubes Thawing:** Store in the fridge ≤40°F (4°C) for 90 minutes before use, and they are ready to use immediately. After freezing at 0°F (-18°C) for 45 minutes, thaw at 59°F-77°F (15°C-25°C) for 15 minutes before use.

TIME OF USE

The following test data is based on laboratory experiments conducted at specific test temperatures, and the data is for reference only regarding human use.

59°F/15°C PHASE CHANGE MATERIAL

Model	Number of ice bars	Ambient temperature	(35.6°F - 77°F) 2°C-25°C insulation duration
7.78 OZ	1	104°F/40°C	205mins/3.4h
	2		360mins/6.0h
	3		474mins/7.9h
	4		504mins/8.4h
7.78 OZ	1	95°F/35°C	275mins/4.6h
	2		496mins/8.3h
	3		604mins/10.1h
	4		736mins/12.3h
7.78 OZ	1	86°F/30°C	368mins/6.1h
	2		603mins/10.1h
	3		914mins/15.2h
	4		1124mins/18.7h

WORKING PRINCIPLE

In a freezer (0°F/-18°C), the biogel ice tubes solidifies, storing cooling energy. Inside the insulin cooler, vacuum insulation and low-conductivity materials maintain a stable temperature, while a reflective coating minimizes heat loss, ensuring insulin stays properly stored.

USAGE SCENARIOS

The insulin cooler cup is ideal for travel, outdoor activities, and daily use, keeping medication at a stable temperature to ensure its effectiveness.

CAUTIONS

- **Clean Regularly:** Use warm water and mild detergent; check and replace seals as needed.

- **Handle Carefully:** Avoid sunlight, high heat, and impacts to the liner.
- **Inspect biogel ice tubes :** Check for damage or leaks before use; replace if needed.
- **Prevent Frostbite:** Wrap ice packs in a cloth during cold therapy.
- **Use Efficiently:** Adjust the number of ice packs based on cooling needs.
- **Replace Periodically:** Swap out biogel ice tubes as their cooling efficiency decreases.

WARNINGS

Discard Spoiled Insulin: Stop using insulin if it looks cloudy, discolored, or abnormal, and consult a doctor.

Keep Away from Children: Store the cooler out of children's reach to prevent accidents.

Check for Allergies: If allergic to materials like stainless steel or rubber, consult a doctor or use alternatives.

CONTRAINDICATIONS

- **No Microwave Heating:** Do not heat the cooler or ice packs in a microwave.
- **No Chemical Cleaners:** Avoid using chemical solvents for cleaning.
- **No High-Temperature Sterilization:** Do not sterilize the cooler or ice packs at high temperatures.
- **Do Not Ingest:** The cooling liquid inside biogel ice tubes contains chemicals. If ingested, seek medical attention immediately.
- **Avoid Electrical Appliances:** Melting biogel ice tubes may release water that could cause electrical short circuits. Keep them away from electronics.
- **No Use with Flammable Items:** Melting biogel ice tubes may release water vapor, posing risks with flammable or explosive items. Avoid using them in such cases.

FAQS

Q: Can the insulin cooler cup store other medications?

A: Yes, medications requiring a storage temperature of 35.6°F-46.4°F (2°C-8 °C) can be stored temporarily. However, long-term storage is not recommended as different medications may have specific storage requirements. Always follow the medication's storage guidelines.

Q: Can I clean the cooler cup with boiling water?

A: It's not recommended to use boiling water for cleaning as high temperatures may affect the cooler cup's insulation performance and materials. Use warm water and mild detergent instead.

Q: What factors affect the insulation duration?

A: The insulation time depends on factors such as external temperature, the initial temperature of the cooler cup, and how frequently the lid is opened. Higher external temperatures and frequent lid openings may reduce insulation time.

PRODUCT INFORMATION

Brand: TEMPK™

Name: Insulin Cooler

Size: 46*46*208mm

Liner material: SUS304


Plastic parts: Food-grade PP

Sealing ring: Food-grade silicone

Producer: Shanghai Huizhou Industrial Co., Ltd.

Origin: Made in China

PACKAGE CONTAINS

x1

7.78 OZ Insulin Cooler

x5

Biogel ice tubes

CONTACT US

If you need assistance, please contact us.

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59°F/15°C PHASE CHANGE MATERIAL			
Model	Number of ice tubes	Ambient temperature	(35.6°F - 77°F) 2°C-25°C insulation duration
15.9 OZ	2	95°F/35°C	372mins/6.2h
	4		606mins/10.1h
	6		672mins/11.2h
	8		816mins/13.6h
15.9 OZ	2	86°F/30°C	510mins/8.5h
	4		876mins/14.6h
	6		> 1200mins/ > 20h
	8		> 1200mins/ > 20h

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CAUTIONS

- **Clean Regularly:** Use warm water and mild detergent; check and replace seals as needed.
- **Handle Carefully:** Avoid sunlight, high heat, and impacts to the liner.
- **Inspect biogel ice tubes:** Check for damage or leaks before use; replace if needed.
- **Prevent Frostbite:** Wrap ice packs in a cloth during cold therapy.
- **Use Efficiently:** Adjust the number of ice packs based on cooling needs.
- **Replace Periodically:** Swap out biogel ice tubes as their cooling efficiency decreases.

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FAQS


- Q:** Can the insulin cooler cup store other medications?
- A:** Yes, medications requiring a storage temperature of 35.6°F-46.4°F (2°C-8°C) can be stored temporarily. However, long-term storage is not recommended as different medications may have specific storage requirements. Always follow the medication's storage guidelines.
- Q:** Can I clean the cooler cup with boiling water?
- A:** It's not recommended to use boiling water for cleaning as high temperatures may affect the cooler cup's insulation performance and materials. Use warm water and mild detergent instead.
- Q:** What factors affect the insulation duration?
- A:** The insulation time depends on factors such as external temperature, the initial temperature of the cooler cup, and how frequently the lid is opened. Higher external temperatures and frequent lid openings may reduce insulation time.

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Fabricante: Shanghai Huizhou Industrial Co., Ltd.